



Calibration solutions for milk cryoscopes

Designation	Freezing point °C
Cryoscopy solution 00,0	0,000 ± 0,002°C
Cryoscopy solution 40,8	- 0,408 ± 0,002°C
Cryoscopy solution 51,2	- 0,512 ± 0,002°C
Cryoscopy solution 52,0	- 0,520 ± 0,002°C
Cryoscopy solution 55,7	- 0,557 ± 0,002°C
Cryoscopy solution 60,0	- 0,600 ± 0,002°C
Cryoscopy solution 80,0	- 0,800 ± 0,003°C